

GEBRAUCHSANWEISUNG

FROWEIN GMBH & CO. KG

SCHWABEX-PRESS




BIOZIDE



Lightweight and easy to handle gel-applicator to apply SchwabEX-gel, SchwabEX-prime, JuvenEX-gel, Monitoring-gel and other gel formulations.



Detailed view with gel-cartridge.



Slide-in the press-out piston (1) from the front to the pistol grip (2), thereto push down the unlocking device (3) and push back the entire press-out piston. Insert cartridge and lock it in 90-degree rotation in cartridge holder (4).